

Exclusion diagnostic

D-Dimer assay

STA-Liatest D-Di Plus

Etiologic diagnostic: screening tests

Antithrombin
assay

STA-Stachrom AT III

Protein C
assay

STA-Staclot
Protein C

or

STA-Stachrom
Protein C

Protein S
assay

STA-Staclot
Protein S

+

STA-Stachrom
Protein C

ou

Asserachrom
Free Protein S

Activated Protein C
resistance assay

STA-Staclot APC-R

Etiologic diagnostic: additional tests

Protein C
assay

Asserachrom
Protein C

Protein S
assay

Asserachrom
Total Protein S